V-CHEK™ CE

IVD

V-CHEK™ SARS-CoV-2 Rapid Antigen Test Kit

For Professional Use / Catalogue No : SCVCO2

Figure 1. Sampling


Figure 2. Processing the sample


Figure 3. Interpretation of test results

Figure 4. Line Interpretation

COMPONENTS	1 Test / box	25 Test / package
Tester (for single use)	1 test cassette (1 test / bag x 1 bag)	25 Test cassette (1 test / bag x 25 bags)
Tampon (for single use)	1 extraction tube, 500 µL extraction tampon	25 extraction tubes, each of which is 500 µL extraction tampon
Sampling sticks (for single use)	1 sterile, disposable sampling swab	25 sterile, disposable sampling stick
Package supplement	1 instruction sheet	1 instruction sheet

Table 1. Materials and Components

SARS-CoV-2 Rapid Antigen Test Kit	RT-Test result in PCR comparison		
	Positive	Negative	Total
Positive	674	5	679
Negative	18	572	590
Total	692	577	1269

Sensitivity: %97.40: (674 / 692), (%95 GA: 95.92 - 98.45)
Specificity: %99.13: (572 / 577), (%95 GA: 97.99 - 99.72)
Accuracy: %98.19: (674 + 572) / (1269)

Table 2. Performance Data (Nasopharyngeal swab)

SARS-CoV-2 Rapid Antigen Test Kit	RT-Test result in PCR comparison		
	Positive	Negative	Total
Positive	156	0	156
Negative	7	110	117
Total	163	110	273

Sensitivity: %95.71: (156 / 163), (%95 GA: 91.35 - 98.26)
Specificity: %100: (110 / 110), (%95 GA: 96.70 - %100.00)
Accuracy: %97.44 (156 + 110) / 273

Table 3. Performance Data (Oropharyngeal swab)

SARS-CoV-2 Rapid Antigen Test Kit	RT-Test result in PCR comparison		
	Positive	Negative	Total
Positive	580	5	585
Negative	20	495	515
Total	600	500	1100

Sensitivity: %96.66: (580 / 600), (%95 GA: 94.90 - 97.95)
Specificity: %99.00: (495 / 500), (%95 GA: 97.68 - 99.67)
Accuracy: %97.72: (580 + 495) / (1100)

Table 4. Performance Data (Nasal swab)

PURPOSE OF USAGE

The V-CHEK™ SARS-CoV-2 Rapid Antigen Test Kit is a lateral flow swab test designed for the qualitative/non-quantitative in vitro detection of SARS-CoV-2 nucleocapsid antigen in nasal, nasopharyngeal and oropharyngeal swab samples.

This test is only suitable for use in clinical laboratories or in close-to-patient examinations by professional users and provides assistance in diagnosing of a SARS-CoV-2 infection. The test is not designed for personal use. The test should not be used as the sole criterion to diagnose the pneumonia arising from SARS-CoV-2 infection and to reach the conclusion that there is any pneumonia. A negative test result does not indicate the absence of SARS-CoV-2 infection. It is recommended to combine the patient's clinical findings and other laboratory tests for a comprehensive analysis of the disease. The test is suitable for patients of all ages.

SUMMARY AND EXPLANATION
The novel coronavirus SARS-CoV-2 is a positive-stranded RNA virus and is a genus of β-Coronaviruses. COVID-19 is an infectious disease of the acute respiratory tract to which humans are susceptible. Currently, patients infected by SARS-CoV-2 are the main source of infection; Asymptomatic infected persons can also transmit the virus. According to current epidemiological research, the incubation period varies between 1 -14 days, mostly 3 - 7 days. The main symptoms are fever, fatigue, loss of smell and dry cough. A small number of nasal congestion, runny nose, sore throat, muscle pain and diarrhea are seen.

PRINCIPLE OF TEST
This reagent uses a double swab method with antibody for the qualitative detection of SARS-CoV-2 nucleocapsid antigen. During the application of test, a monoclonal anti-SARS-CoV-2 antibody labelled with colloidal gold are connected to the SARS-CoV-2 antigen in the sample. This reaction complex proceeds chromatographically on the nitrocellulose membrane and is connected to the pre-coated monoclonal anti-SARS-CoV-2 antibody at the detection zone (T) on the test membrane and creates a red reaction line there. If the sample does not contain SARS-CoV-2 antigen, a red colour reaction line cannot occur in the T zone.

At the same time, a Chicken-IgY-Gold conjugate moves across the membrane during testing and are connected to a pre-coated monoclonal Anti-Chicken-IgY-Antibody in quality control area C, creating a red reaction line there. Regardless of whether the sample to be examined contains the SARS-CoV-2 antigen or not, a red reaction line always occurs in the quality control area (C).

MATERIALS AND COMPONENTS
Materials supplied together with test kits (Page 1, Table 1)
Note: Components of different kit lots should not be mix.
Additional required materials
Personal protective equipment, timers.
Active components of the test Reagents
• mAb anti-COVID-19 antibodies grown in mice
• mAb Anti-Chicken-IgY grown in mice
• mAb Anti-COVID-19 gold conjugated antibody grown in mice
• Purified Chicken-IgY-Gold Conjugate
• Recombinant COVID-19 nucleocapsid protein

STORAGE AND DURABILITY
1. Store the test kit at 2 °C - 30 °C. Do not store or freeze the kit below 2 °C. All components must be brought to room temperature before performing the test.
2. The test cassette should be used within 15 minutes after it is removed from the foil pouch.
3. The kit should not be used after its expiration date. The expiration date is stated on the label/packaging.

TEST PROCEDURE
Before testing, carefully read the instructions for use and follow the below mentioned instructions as described. Make sure that test components are at room temperature when they are used. The testing procedure includes the following steps: sampling, processing the sample, and performing the test.
Caution: Samples should be used as soon as possible after they are collected. After half an hour, the sample should not be used and a new sample should be taken by a new sample stick.
Caution: Samples should not be inactivated.
Caution: The way the sample is taken shows difference between each swab samples. Please perform only one of the specified swab samples (1a - 1c).
1a. Nasopharyngeal spreading: Ask the patient to tilt his/her head back slightly. Then insert slowly the sterile swab with its head first into the nasopharynx until you feel a slight resistance. Turn the swab 3 times in a place near the inner wall of the nasal cavity and carefully remove the swab from the nose. Avoid contact with nasal mucosa when entering and exiting.

or 1b. **Oropharyngeal spreading:** Direct the sterile swab from the side of the uvula towards the back of the pharynx. Rub and twist the swab 10 times along the back of the pharynx and both tonsils. Then remove the stick. Avoid contact of the head of the swab with the tongue during sample collection.
or 1c. **Anterior nasal swab:** Insert the sterile swab into the anterior nasal part and turn the swab 3 times along the inner wall of the nasal cavity. Then remove the stick. (Page 1 Figure 1)
2. Open the cap (drip cap) of the extraction tube and place the used swab into the extraction tube with the swab head first.
3. Turn the swab in the extraction tampon 10 times

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1. Detection limit
At a virus culture concentration of 100 TCID₅₀/mL and more positive value is equal to 95% or greater than it. The detection limit of the SARS-CoV-2 Antigen Test is 100 TCID₅₀/ml
1. Cross-reactivity
The cross-reactivity of the kit was evaluated. The results did not show any cross-reactivity with the following agents at the given concentrations. (Page 13 Table 5)
4. Impeding agents
Test results are not affected by the following agents at the specified concentration:
(Page 14 Table 6)
5. Sensitivity
1. 10 copies of negative and positive samples were tested by using reference materials. The compatibility of negative and positive results is 100%.
2. Three different lots were tested with positive and negative reference materials. The compatibility of negative and positive results is 100%.
6. Hook-Effect
1. No Hook-Effect was found at 0x10 concentration of TCID₅₀/ml SARS-CoV-2.

2. Positive result: If both C and T lines appear after 15-20 minutes, the test result is positive and valid. If your test result is positive, please contact. Contact your family doctor immediately by running the RT-PCR test to confirm the result. To reduce the risk of infection, it is strongly recommended to isolate you and your close contacts quick and to follow standard operating procedures in accordance with current national guidelines and protocols and to receive medical assistance.
2. Negative result: If only the C-line becomes visible, but the T-line does not become visible after 15-20 minutes, the test result is negative and valid. If you develop the symptoms of Covid-19, you and your home should be isolated and an RT-PCR test should be done to confirm the result. have an RT-PCR test to confirm the result. You should follow Standard Operating Instructions according to protocol and continue to follow national and local rules and guidelines, including regular hand washing, social distancing and face protection, and see a doctor if necessary.
Invalid result: If no C- line becomes visible after 15-20 minutes, the test result is invalid. If the T-line becomes visible but there is no C-line, the test result is also invalid. In both cases the test should be performed with a new test cassette.

(Page 1 Figure 3, 4)
LIMITATIONS
1. The result of the product should not be considered as a confirmed diagnosis. Evaluation of test results should be done together with RT-PCR results and in addition clinical signs, epidemiological information and more clinical data should be collected.
2. The contents of the kit are designed for the qualitative detection of SARS-CoV-2 antigens from nasal, oropharyngeal and nasopharyngeal swabs. Other sample types cannot be used.
3. Test performance is dependent on the amount of virus (antigen) in the sample and it may or may not correlate with the results of a virus culture performed on the same sample.
4. A negative test result may occur if the antigen concentration in a sample is below the detection limit of the test or if the sample has been improperly collected or transported.
5. If test procedure is followed, this may adversely affect test performance and/or invalidate the test result.
6. If the result is read before 15 minutes, this can lead to a false negative result. If the result is read after 20 minutes, this may cause a false positive result.
7. A positive test result does not exclude the possibility of infection together with other pathogens.
8. Negative test results do not eliminate other viral or bacterial infections that are not SARS.
9. Negative results should be considered as assumption and confirmed by a molecular test.
10. A positive result should not be taken as a confirmed diagnosis. Evaluation should be done in conjunction with clinical symptoms and other diagnostic methods.
11. False negative/positive results may occur if a sample is improperly collected, transported or processed.
12. There is also a possibility that the test will give a positive result even if there is no covid infection (false positive)
13. Improper sample collection, use of other reagents that do not match the test kits, and incorrect application of test protocols for the elution process may result in false positive/negative results.

PERFORMANCE DATA
1. Clinical validation
The clinical performance of V-CHEK has been determined by comparing the V-CHEK™ SARS-CoV-2 Rapid Antigen Test with a nasopharyngeal RT-PCR test as nasopharyngeal. Samples are taken 7 days after symptoms appear.
a) Nasopharyngeal swab
The performance of the SARS-CoV-2 Rapid Antigen Test Kit was evaluated by means of using 1269 nasopharyngeal swabs received from patients. (Page 1 Table 2)
b) Oropharyngeal swab
The performance of the SARS-CoV-2 Rapid Antigen Test Kit was evaluated by means of using 273 oropharyngeal swabs received from patients. (Page 1 Table 3)
c) Nasal swab
The performance of the SARS-CoV-2 Rapid Antigen Test Kit was evaluated using 1100 nasal specimen (nasal) swabs from patients. (Page 1 Table 4)

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